



## **Trackmaster Treadmill mod. TMX425**

### **FEATURES**

We built the Trackmaster TMX425 to be the best stress treadmill in the industry.

We did so by making it smoother, quieter and more reliable than any competitor, while centering the design on the absolute safety and comfort of the patient.

The treadmill is self-calibrating, virtually maintenance-free and the safety features include a low profile deck, incredibly smooth belt operation and an easy-to-reach emergency stop switch.

The large 22" x 63" deck absorbs impact while providing plenty of room to run or walk, and ergonomically designed handrails inspire patient confidence.

Available manual and fully programmable rail-mounted controllers are simple and intuitive, allowing the treadmill to be used independently from the stress system.

***Automatic belt alignment***

***Smooth & Quiet***

***Safe & comfortable***

***Rugged & reliable***

## SPECIFICATIONS

## TMX425

Patented V-Belt Track provides smooth, safe and comfortable operation  
 Whisper quiet operation  
 Large running surface  
 Low profile deck 7 1/2" from floor  
 Emergency stop button control



182kg (400lbs) weight capacity  
 Double-sided running deck absorbs and dampens the shock while providing extended usage  
 Fully integrated Integral RS232 interface  
 Seamless communication protocol provided  
 Optional paediatric hand rail  
 Rugged & reliable  
 All metal chassis  
 Simplified drive assembly provides easy access for maintenance  
 Self-adjusting automatic belt alignment

**Specifications****Treadmill**

182kg (400lbs) capacity

**Drive System**

Standing start  
 Self-aligning running belt  
 Heavy duty 2.2 HP AC motor  
 110V, 60 Hz, 20A  
 208V/220V, 60 Hz, 15A  
 Dedicated circuit for each treadmill

**Speed Range**

0.8 – 19.3 kph (0.5 – 12 mph) for 220V motor  
 0.8 - 16.1 kph (0.5 - 10 mph) for 110V motor  
 Auto Speed Reset (ASR)

**Elevation range**

0-25% (0-22.5°)

**Floor surface requirements**

0.82x2.24m (33"x88") footprint

**Running surface**

Shock absorbing  
 0.56x1.60m (22"x63")  
 0.19m (7 1/2") from floor  
 Double-sided running deck

**Safety features**

Emergency stop button  
 Safety handrails

**Shipping weight**

227kg (500lbs) uncrated on pallet  
 273kg (600lbs) crated